CE 0123

DE Tageslichtlampe Gebrauchsanweisung.....	2	Electromagnetic Compatibility Information	65
EN Daylight therapy lamp Instructions for use.....	10		
FR Lampe de lumineothérapie Mode d'emploi	18		
S Dagsljuslampa Bruksanvisning	26		
FI Kirkasvalolamppu Käyttöohje	34		
DA Dagslyslampe Betjeningsvejledning	42		
NO Dagslyslampe Bruksanvisning	50		
RU Прибор дневного света Инструкция по применению	58		

Inhalt

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Content

1. About the daylight therapy lamp.....	10	7. Cleaning and care of the unit.....	14
2. Signs and symbols	11	8. Storage	15
3. Notes	11	9. Disposal	15
4. Unit description	13	10. What if there are problems?.....	15
5. Initial use.....	13	11. Technical specifications.....	15
6. Operation	13	12. Warranty/service	16

1. About the daylight therapy lamp

Dear valued customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air.

With kind regards

Your Beurer team

Why use a daylight therapy lamp?

When the hours of daylight are noticeably shorter in autumn and people increasingly stay inside in winter, the effects of a lack of light may become apparent. This is often described as “winter depression”. The symptoms can present themselves in a number of ways:

- Imbalance
- Low mood
- Lack of energy and listlessness
- Generally feeling under the weather
- Increased need for sleep
- Loss of appetite
- Difficulty concentrating

The cause of these symptoms is the fact that light, particularly sunlight, is essential for life and has a direct effect upon the human body. Sunlight indirectly controls the production of melatonin, which is only passed to the blood in darkness. This hormone shows your body that it is time to sleep. Therefore, in months with less sunshine more melatonin is produced, making it difficult to get up in the mornings because your body functions are powered down. Use the daylight therapy lamp immediately after waking up (i.e. as early as possible) to end the production of melatonin and to brighten your mood.

Lack of light also prevents the production of the happy hormone serotonin, which is said to significantly influence our well-being. The application of light thus yields quantitative changes to hormones and neurotransmitters in the brain that have an effect on our activity levels, our feelings and our well-being. To compensate for such a hormonal imbalance, daylight therapy lamps can create a suitable replacement for natural sunlight.


In the medical field, daylight therapy lamps are used as compensation for the effects of a lack of light. Daylight therapy lamps simulate daylight over 10,000 lux. This light can influence the human body and be used as a treatment or as a preventative measure. Normal electric light, however, is not sufficient to influence the hormonal balance. This is because in a well-lit office, the light intensity is just 500 lux, for example.

Included in delivery:

- Daylight therapy lamp
- These operating instructions

2. Signs and symbols

The following symbols appear in these instructions or on the device.

	Read the instructions of use		Permissible operating temperature and humidity
 Warning	Warning instruction indicating a risk of injury or damage to health	 Note	Note on important information
 Important	Safety note indicating possible damage to the unit/accessory		Manufacturer
	Serial number		Device in protection class 2
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products		OFF / ON
	Protected against solid foreign objects 12.5 mm in diameter and larger, and against vertically falling drops of water		Dispose of packaging in an environmentally friendly manner
 Storage	Permissible storage temperature and humidity		

3. Notes

Please read these instructions for use carefully and keep them for later use. Be sure to make them accessible to other users and observe the information they contain.

Safety notes

Warning

- The daylight therapy lamp is only intended for radiation on human bodies.
- The device is not intended for commercial or clinical use; it is designed exclusively for personal use at home.
- The device is suitable for health applications at home. No special training is required to use the device.
- Before use, ensure that there is no visible damage to the device or accessories and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Service address.
- Ensure that the daylight therapy lamp is positioned on a stable base.
- The device must only be connected to the mains voltage that is specified on the type plate.
- Never submerge the device in water or other liquids and do not use it in the bathroom.
- It must not be used on people with disabilities, children under 3 years of age or people with reduced sensitivity to heat (e.g. people with skin alternations due to illness), after taking pain relief medication, alcohol or drugs.
- Keep packaging material away from children (risk of suffocation!).
- Do not allow children to play with the device.
- If the device is warm, do not cover it or pack and store it.
- Always unplug the mains part and allow the device to cool down before touching it.
- Do not touch the device with wet hands when it is plugged in and do not allow water to spray on the device. Only operate the device if it is completely dry.
- Ensure that you only insert and remove the mains part with dry hands and that you only press the ON/OFF button with dry hands.
- Keep the mains cable away from hot objects and naked flames.
- Protect the device from heavy impact.
- Do not pull the mains part out of the socket using the mains cable. Hold the mains part directly when pulling the mains part from the socket.

- Do not use the device if it shows signs of damage or does not function properly. In such cases, contact Customer Services.
- If the mains connection cable of this device is damaged, it must be disposed of. If it cannot be removed, the device must be disposed of.
- The plug is used to disconnect the power supply. Position the device so that you can easily unplug the plug from the power supply at all times.
- Disconnection from the power supply network is only guaranteed when the mains part is unplugged.
- Do not use the device in the presence of flammable anaesthetic gas connections with air, oxygen or nitro-gen oxide.
- No calibration and no preventive checks or maintenance need to be carried out on this device.
- You cannot repair the device. The device contains no parts that you can repair.
- Do not make any changes to the device without the manufacturer's permission.
- If the device has been changed, thorough tests and checks must be carried out to ensure the continued safety of further use of the device.

General notes

Important

The following conditions may be contraindication for use of this appliance:

- suffer(ed) from severe depression;
- recent eye surgery or a diagnosed eye condition whereby doctor has advised you to avoid bright light;
- skin sensitive to light, such as systemic lupus erythematosus.
- Always consult a doctor before using the daylight therapy lamp if you are taking medication such as pain relief medication, medication to reduce high blood pressure or antidepressant medication.
- Diabetics and people who suffer from retinal diseases must be examined by an optician before using the daylight therapy lamp.
- Please do not use the device if you suffer from an eye disease such as cataracts, glaucoma, diseases of the optic nerve or inflammation of the vitreous body.
- If you have health concerns of any kind, consult your GP!
- Always consult a doctor before using the daylight therapy lamp if you have a strong sensitivity to light or you are prone to migraine attacks.
- Remove all packaging material before using the device.
- Light sources are excluded from the warranty.

EMC information

Warning

- Use of this appliance adjacent to or stacked with other device should be avoided because it could result in improper operation. If such use is necessary, this appliance and the other device should be observed to verify that they are operating normally. If the appliance does not light on or is flickering, it means that the appliance is not operating normally.
- Do not change any cable without authorization of the manufacturer. Use of cables other than those specified or provided by the manufacturer of this appliance could result in increased electromagnetic emissions or decreased electromagnetic immunity of this appliance and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm to any part of this appliance, including power cable. Otherwise, degradation of the performance of this equipment could result.
- If the appliance does not light on or is flickering due to the electromagnetic disturbances, do not use it and try to place the appliance in other location.
- For maintaining the basic safety and essential performance with regard to electromagnetic disturbances for the expected service life, please take precaution if the use location is near (e.g. less than 1,5 km from) from AM, FM or TV broadcast antennas.

Repairs

Important


- Do not open the device. Do not attempt to repair the device yourself. This can result in serious injury. Failure to comply with this instruction will void the warranty.
- For repairs, please contact Customer Services or an authorised retailer.

- Do not use the device if it is damaged. If the device is not working, try to solve the problem using the information in chapter “What if there are problems?”. If you are unable to solve the problem, contact Customer Services or an authorised retailer.

4. Unit description

Overview

1. Fluorescent screen
2. Rear of the device housing
3. Support feet
4. On/off button
5. Mains cable with plug
6. Stand
7. Hole for wall-mounting fixture

5. Initial use

Take the device out of the plastic wrapping. Check the device for damage or faults. If you notice any damage or faults on the device, do not use it and contact customer services or your supplier.

Positioning the device

Place the device on an even surface. The position should be chosen to ensure a distance of between 20 cm and 45 cm between the user and the device. The lamp is most effective at this distance.

Mounting the device on the wall

Alternatively, you can also mount the device on the wall. This requires a solid wall with a stable suspension device (for concrete walls e.g. a 25 mm long M4 stainless steel self-tapping screw and wall plug). Only mount the device on the suspension device at the hole (located centrally on the rear of the device). Self-tapping screw and wall plug are not included in delivery.

Important:

If the area around the hole in the rear of the device is cracked or damaged, do not mount the device on the wall! For help, please contact customer services or your retailer.

Note:

- Ensure that there is a mains socket close to the set-up area.
- Ensure that the mains part plug is completely inserted in the mains socket.
- Arrange the mains cable so that no-one will trip over it.

6. Operation

1	Switching on the lamp ▶ Press the On/Off button.
----------	--

2	<p>Enjoying the light</p> <p>Sit as close as possible to the lamp, between 20 cm and 45 cm. You can go about your normal activities while using it. You can read, write, make telephone calls, etc.</p> <ul style="list-style-type: none"> • Every so often, briefly look directly into the light, since the effect results from the eyes/ retina receiving the light. • Use the daylight therapy lamp as often as you want. However, the treatment is at its most effective if you carry out the light therapy according to the prescribed time period of at least 7 successive days. • The most effective time of day for the treatment is between 6am and 8pm and we recommend that it is used for 2 hours per day. • Do not look directly into the light the whole time you are using it, however, since that could lead to overstimulation of the retina. • Start with shorter periods of illumination, which you increase over the course of a week. • Do not leave the daylight therapy lamp unsupervised if it is being used. <p>Note:</p> <p>After the first few applications you may experience side effects such as pain in the eyes and head or nausea. This should disappear in further sessions as the nervous system becomes accustomed to the new stimuli. Refrain from using the daylight therapy lamp for a few days until the side effects have subsided, then try to use the device again.</p>												
3	<p>Important instructions</p> <p>When using the lamp, maintain the recommended distance of 20 - 45 cm between your face and the lamp. The duration of the application depends on the distance:</p> <table border="1" data-bbox="143 628 771 775"> <thead> <tr> <th>Lux</th> <th>Distance</th> <th>Application duration</th> </tr> </thead> <tbody> <tr> <td>10.000</td> <td>approx. 20 cm</td> <td>0,5 hour</td> </tr> <tr> <td>5.000</td> <td>approx. 30 cm</td> <td>1 hour</td> </tr> <tr> <td>2.500</td> <td>approx. 45 cm</td> <td>2 hours</td> </tr> </tbody> </table> <p>Basically, the closer you are to the source of light, the less time you should use it.</p>	Lux	Distance	Application duration	10.000	approx. 20 cm	0,5 hour	5.000	approx. 30 cm	1 hour	2.500	approx. 45 cm	2 hours
Lux	Distance	Application duration											
10.000	approx. 20 cm	0,5 hour											
5.000	approx. 30 cm	1 hour											
2.500	approx. 45 cm	2 hours											
4	<p>Enjoying light over longer periods</p> <p>Repeat your use of the lamp in the dark seasons for at least 7 successive days, or even longer, depending on your individual needs. If possible, conduct the treatment in the morning hours.</p>												
5	<p>Switching off the lamp</p> <p>► Press the ON/OFF button. The LEDs switch off. Unplug the mains part from the mains socket.</p> <p>CAUTION!</p> <p>The lamp remains warm after use. Allow the lamp to cool down sufficiently before putting it away and/or packing it away warm.</p>												

7. Cleaning and care of the unit

The device should be cleaned from time to time.

Important

- Ensure that no water gets inside the device.
The device must be switched off, disconnected from the mains and allowed to cool down each time before cleaning.
- Do not clean the device in a dishwasher.
Use a slightly damp cloth to clean the device.
- Do not use any abrasive cleaning products and never submerge the device in water.
Do not touch the device with wet hands when it is plugged in and do not allow water to spray on the device. Only operate the device if it is completely dry.

8. Storage

If you are not going to use the device for an extended period of time, disconnect it and store it in a dry place, out of the reach of children.

Follow the storage instructions provided in the “Technical data” chapter.

9. Disposal

For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

10. What if there are problems?

Problem	Possible cause	Solution
Device does not light up	On/off button switched off	Switch on the on/off button.
	No electricity	Connect the mains plug correctly.
	No electricity	The mains cable is defective. Contact customer services or your retailer.
	LEDs have reached the end of their service life. LEDs faulty.	For repairs, please contact Customer Services or an authorised retailer.

11. Technical specifications

Type	MDLEDV2
Model	TL41
Dimensions (WHD)	270 x 207 x 80 mm
Cable length	max. 2.000 mm
Weight	approx. 645 g
Light	LEDs
Illuminance	10,000 lux (at a distance of approx. 20 cm)
Mains connection	220-240 V ~, 50 Hz, 0.3 A
Operating conditions	0°C to +40°C, 30-90% relative humidity, 700-1060 hPa ambient pressure
Storage conditions	-10°C to +60°C, 30-90% relative humidity, 700-1060 hPa ambient pressure
Transport conditions	-10°C to +60°C, 30-90% relative humidity, 700-1060 hPa ambient pressure
Product classification	Protection class II, IP21, no AP or APG, continuous operation
Location of emission aperture	Screen (around approx. 440 cm ²)
Color temperature of LEDs	6.500 Kelvin
LED lifetime	approx. 10,000 hours
Risk group of light source equipment	risk group 1 classified by IEC 60601-2-57: 2011
Product lifetime	5 years

Maximum output of optical radiation for the TL 41

Radiance Output	Risk group classified by IEC 60601-2-57: 2011	Max. value
E _{UVA} : Eye UV-A	Exempt Group	$1.6 \times 10^{-4} \text{ W}\cdot\text{m}^{-2}$
ES: Actinic UV skin & eye	Exempt Group	$1.5 \times 10^{-4} \text{ W}\cdot\text{m}^{-2}$
EIR: Infrared radiation hazard exposure limits for the eye	Exempt Group	$5.7 \times 10^{-2} \text{ W}\cdot\text{m}^{-2}$
EH: Thermal hazard for the skin	Exempt Group	$3.1 \times 10^1 \text{ W}\cdot\text{m}^{-2}$
LB: Blue light	Risk Group 1	$1.1 \times 10^3 \text{ W}\cdot\text{m}^{-2}\cdot\text{sr}^{-1}$
LR: Retinal thermal	Exempt Group	$2.4 \times 10^2 \text{ W}\cdot\text{m}^{-2}\cdot\text{sr}^{-1}$

Subject to technical modifications.

The serial number is located on the device or in the battery compartment.

This device meets the requirements of the European Directive 93/42/EEC and 2007/47/EC Medical Devices as well as the Medical Devices Act.

This device complies with European Standard EN 60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility. Please note that portable and mobile HF communication systems may interfere with this unit. More details can be requested from the stated customer service address or found at the end of the instructions for use.

12. Warranty/service

Beurer GmbH, Söflinger Straße 218, 89077 Ulm, Germany (hereinafter referred to as “Beurer”) provides a warranty for this product, subject to the requirements below and to the extent described as follows.

The warranty conditions below shall not affect the seller’s statutory warranty obligations which ensue from the sales agreement with the buyer.

The warranty shall apply without prejudice to any mandatory statutory provisions on liability.

Beurer guarantees the perfect functionality and completeness of this product.

The worldwide warranty period is 3 years, commencing from the purchase of the new, unused product from the seller.

The warranty only applies to products purchased by the buyer as a consumer and used exclusively for personal purposes in the context of domestic use.

German law shall apply.

During the warranty period, should this product prove to be incomplete or defective in functionality in accordance with the following provisions, Beurer shall carry out a repair or a replacement delivery free of charge, in accordance with these warranty conditions.

If the buyer wishes to make a warranty claim, they should approach their local retailer in the first instance: see the attached “International Service” list of service addresses.

The buyer will then receive further information about the processing of the warranty claim, e.g. where they can send the product and what documentation is required.

A warranty claim shall only be considered if the buyer can provide Beurer, or an authorised Beurer partner, with

- a copy of the invoice/purchase receipt, and
- the original product.

The following are explicitly excluded from this warranty:

- deterioration due to normal use or consumption of the product;
- accessories supplied with this product which are worn out or used up through proper use (e.g. batteries, rechargeable batteries, cuffs, seals, electrodes, light sources, attachments and nebuliser accessories);
- products that are used, cleaned, stored or maintained improperly and/or contrary to the provisions of the instructions for use, as well as products that have been opened, repaired or modified by the buyer or by a service centre not authorised by Beurer;
- damage that arises during transport between manufacturer and customer, or between service centre and customer;
- products purchased as seconds or as used goods;
- consequential damage arising from a fault in this product (however, in this case, claims may exist arising from product liability or other compulsory statutory liability provisions).

Repairs or an exchange in full do not extend the warranty period under any circumstances.

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ELECTROMAGNETIC COMPATIBILITY INFORMATION

The "TL 41" has been designed and manufactured to conform to the highest standards of safety and performance including electromagnetic compatibility (EMC).

The appliance is complied with the following emission and immunity standard, and test level:

Emissions Test Standard	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Comply

Guidance and manufacturer's declaration – electromagnetic immunity							
Radiated RF IEC61000-4-3 (Test specification for ENCLOSED PORT IMMUNITY to RF wireless communications equipment	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380 - 390	TETRA 400	Pulse modulation 18 Hz	1,8	0.3	27
	450	430 - 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704 - 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0.3	9
	745						
	780						
	810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
	870						
	930						
	1720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0.3	9
	5240						
	5785						

Immunity Test Standard	IEC 60601-1-2 Test Level	Compliance Level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8kV contact ± 2kV, ± 4kV, ± 8kV, ± 15kV air	± 8 kV contact ± 2kV, ± 4kV, ± 8kV, ± 15kV air
Electrical fast transient/burst IEC 61000-4-4	± 2kV, for power supply lines ± 1kV signal input/ output 100 kHz repetition frequency	±2 kV, for power supply lines Not applicable 100 kHz repetition frequency
Surge IEC 61000-4-5	± 0.5kV, ± 1kV differential lines ± 0.5kV, ± 1kV, ± 1kV common mode	± 0.5kV, ± 1kV differential lines Not Applicable
Voltage dips, interruptions and variations on power supply input lines IEC 61000-4-11	0 % U _T ; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % U _T ; 1 cycle and 70 % U _T ; 25/30 cycles; Single phase: at 0°. 0 % U _T ; 250/300 cycle	0 % U _T ; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % U _T ; 1 cycle and 70 % U _T ; 25/30 cycles; Single phase: at 0°. 0 % U _T ; 250/300 cycle
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz

Note: U_T is the a.c. mains voltage prior to application of the test level.

Zero-Plus International Limited
Rm 1004, 10/F, Join-in Hang Sing Centre, 71-75 Container, Port Road
Kwai Chung, New Territories, Hongkong

Distributed by

BEURER GmbH • Söflinger Str. 218 • 89077 Ulm (Germany)
www.beurer.com • www.beurer-gesundheitsratgeber.com
www.beurer-healthguide.com

Globalmind Consumer Electronics GmbH
Ernst-Mantius-Str.11,
21029 Hamburg, Germany